

Sensors	No.	No.
mini sensor :	1	3
Due sensor (A)	2	6
Due sensor (B)	2	6
big sensor	11	13

Test goal :

1. sensor different

2. wall noise


3. placement error

4. Due for both muscle ??

Test_0 : not attached to skin for reference

Test_1 : Due (A) + big sensor

	FDE	BCP
L	2(A)	11
R	6(A)	13

* EMG reading channel  on Trigon contract
10 rest + 10 MVC

Test_1-inverse: Due (A) + big sensor

	FDE	BCP
L	6(A)	13
R	2(A)	11

10 rest + 10 MVC

Test_2

	Due (B)	mini sensor
	FDE	BCP
L	2(B)	1
R	6(B)	3

10 rest + 10 MLC

Test_2_inverse

	Due (B)	mini sensor
	FDE	BCP
L	6(B)	3
R	2(B)	1

10 rest + 10 MLC

~~Test_3~~ (for placement check)

	mini	Due (A)
L	1	2(A)
R	3	6(A)

10 rest + 10 MLC

Test_4 (for Due test) & same ref point

	Due (B)	Due (A)
L	6(B)	6(A)
R	2(B)	2(A)

10 rest + 10 MLC